

gehealthcare.com




EMPOWERING EXCEPTIONAL CARE

Every pelvic exam brings the opportunity to help ensure the future health of your patients


SIMPLE, SOPHISTICATED IMAGING

2D imaging optimized for clarity and detail with minimal fine-tuning


Obtain views not possible with 2D scanning for a more complete evaluation with **3D Imaging**




Improve contrast resolution and visualization of anatomy in any image plane, even irregularly shaped structures **Advanced VCI with OmniView**



Experience a new standard of color Doppler with **Radiant** flow – Delivering easy, fast visualization of blood flow*



Expand the range of visible blood flow to include low velocities with **Slow** flow HD to visualize blood perfusion in very small vessels*



ENHANCED COMFORT - INNOVATIVE TOOLS

3D Hysterosalpingo-Foam Sonography (HyFoSy) and HyCoSy (Hysterosalpingo Contrast Sonography) – non-ionizing procedure for tubal patency assessment with 3D display of fallopian tubes**

Display more anatomical information in a single image than with standard field of view (FOV) with **Wide-sector** – available on many abdominal and endocavity probes with FOV maximum of 180 degrees Enhance patient comfort and access to angles more difficult to attain with **Beta View** – enables the user to steer the scan plane via crystal movement rather than probe manipulation

Tomographic Ultrasound Imaging (TUI) makes analysis and documentation of dynamic studies easier with a simultaneous view of multiple slices of a volume data set



STREAMLINED WORKFLOW - MAXIMIZE YOUR DAY

Balance schedules while increasing patient satisfaction with **Scan Assistant**. This flexible, customizable exam protocol tool helps increase exam consistency and productivity while documenting for quality assurance purposes

Simplify identification and reporting of uterine abnormalities with **Uterine Classifications** based on ESHRE, ESGE,

and ASRM guidelines***



Ovarian cancer risk estimate based on IOTA (International Ovarian Tumor Analysis) Simple Rules, IOTA LR2 and IOTA ADNEX are IOTA Group Approved on the Voluson Expert Series

Evaluate risk and assign classification of the endometrial lining based on IETA (International Endometrial Tumor Analysis) and IDEA (International Deep Endometriosis Analysis)



© 2018 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, and Voluson are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

August 2018 JB50957XXad *Available only on Voluson E10 and E8

**HyCoSy/HyFoSy not approved in all countries

***ESHRE – European Society of Human Reproduction and Embryology

ESGE – European Society of Gastrointestinal Endoscopy Endoscopy ASRM – American Society of Reproductive Medicine